

Roydent: *Your 1st Choice in Endo*

Endo Forceps



Post & Wedge Forceps – Item R5716

- Ideal for holding posts during placement
- Excellent for placing and removing wedges
- Concave and serrated jaws for firm grip
- Grooved to hold suture needle (similar to Steiglitz Archer 45° forceps)
- Manufactured from German stainless steel



90° Micro Endo Forceps – Item R5606

- Serrated jaws for firm grip
- Ideal for removing silver points
- Straight 8.5cm shaft with 7.0mm beaks
- Manufactured from German stainless steel



45° Advanced ME Forceps – Item R6661

- Ideal for use during endodontic procedures in removing small objects from tight areas
- Advanced design offers a much shorter instrument
- Angled 4.2cm shaft with 8.0mm beaks
- Manufactured from German stainless steel

 **ROYDENT**[™]
DENTAL PRODUCTS

1-800-992-7767
www.roydent.com

Did you know...?


Roydent offers a family of compatible products for **Root Canal Treatment**.

Sodium Hypochlorite (NaOCl) 5 to 5½ % Solution



C.L. — Canal Lubricant Chelates and lubricates canals

Endo-Cleanse (EDTA) 17% solution aids in removing smear layer



Item Number	Description
PL209	EDTA - 8 oz. bottle
PL208	EDTA - 2 oz. bottle
703	EDTA - 28 pipettes
103	C. L. - 20 pipettes
101	C. L. - 150 pipettes
704	NaOCl - 28 pipettes

ROYDENT™
DENTAL PRODUCTS

1-800-992-7767
www.roydent.com