

LOYALTY



LIMITED TIME PROMOTIONAL OFFERS!

*Mix & Match: C-Files • K-Files • Flexicut Files
Hedstroms • Reamers Gates Gliddens Drills
Peeso Reamer Drills

Purchase 25 packs*, get

5 Packs FREE!

Promo Code RIH

Purchase 10 packs*, get

2 Packs FREE!

Promo Code RII


Purchase 5 packs*, get

1 Pack FREE!

Promo Code RIJ

See back for **ANOTHER** special
introductory offer!

Roydent Dental Products prides itself on long standing relationships with dealer partners, universities & students, end user clinicians and staff which also creates loyalty. Customer loyalty has catapulted Roydent into one of the top three hand file companies in the USA. We strive to deliver products and customer service that you know and trust. We are excited to share our Loyalty Purchasing Program with you, please call us for details.



FILE	TIP	DESIGN
C-FILES	Non-Cutting tip	Designed for initial access and breaking through calcification.
FLEXICUT FILES	Non-Cutting tip	Triangular blade for increased flexibility. Designed for full instrumentation.
K-FILES	Cutting tip	Designed for full instrumentation.
NITI FILES	Non-Cutting tip	Conforms to the curvature of the canal during instrumentation.

Also available: Hedstroms, Reamers, Mid-Sized K-Files.

ROYDENTTM
DENTAL PRODUCTS
roydent.com • 1.800.992.7767

discover **Aha!**TM

discover

Aha!™

Make a Great First Impression—Every Time.

Introducing the Aha!™ Family of Impression Materials

We understand how frustrating it can be to THINK you have the perfect impression only to have molds chip, break, or tear in processing. With Roydent's new Aha! family of impression materials, you can feel confident knowing you will get a solid, durable impression that will hold up and stay true. You get a great mold AND a happy patient, the FIRST time. Now that truly IS a great impression!



Aha! VPS X-Light Impression Material NEW SKU: K11892R

- Can be used for two-step or one-step putty-wash impression technique hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation.
- High initial hydrophilicity.
- Working time: 2 min, Working time at 35°C/ 95°F: 1 min, Intraoral set time: 2 min, Total set time: 4 min

Aha! Regular Set Bite Registration Material NEW SKU: K11932R

- Syringeable addition-curing elastomeric bite registration materials
- Used for making accurate occlusal registrations.
- Shore-hardness D 43
- Working time: 30 sec Intraoral set time: 1 min 30 sec Total set time: 2 min

Aha! Fast Set Bite Registration Material NEW SKU: K11961R

- Syringeable addition-curing elastomeric bite registration materials
- Used for making accurate occlusal registrations.
- Can be used for applications which require a hard bite registration silicone.
- Shore-hardness D 43
- Working time: 15 sec, Intraoral set time: 45 sec, Total set time: 1 min

Aha! Putty Fast Set NEW SKU: K11121R

- Standard kneadable putty
- Can be used as two-step putty wash in conjunction with Aha! VPS X-Light
- Can be used as one-step putty wash impression
- Working time: 2 min, Intraoral set time: 2 min, Total set time: 4 min

Aha! VPS Heavy Body Impression Material NEW SKU: K13541R

- Can be used for one or two step impression technique
- Functional impressions
- High initial hydrophilicity.
- Working time: 1 min 20 sec Intraoral set time: 2 min Total set time 3 min, 20 sec

SPECIAL INTRODUCTORY OFFER: Buy 2, Get 1 FREE!

Includes all Aha! registration and VPS materials | code RW1

ROYDENT™
DENTAL PRODUCTS

The boutique brand you've come to love & trust.

roydent.com | 1.800.992.7767 | [f](#) [YouTube](#) [in](#)

Offers good July 1 – September 31, 2017. Please fax invoices for redemption no later than October 31, 2017. Redemption Fax: 888-769-3368 or email to promos@roydent.com, Attn: Coupon Coordinator. Please allow 4-6 weeks for delivery of free products. Package quantities must come from one invoice. Limit 5 redemptions per invoice. May not be combined with any other offer. Not available to contract or formulary customers. Available only in the USA. Please note size of free goods desired.