

2Seal

for a better sealed root canal!

→ hermetic sealing

The goal of a root canal filling is to provide biocompatible and long-term hermetic sealing of the entire root canal system. This prevents penetration of liquids and micro-organisms into the canal from the apical and the coronal regions.



2Seal is a root canal sealing material based on epoxy-amine resin chemistry and consisting of two easy to mix pastes. In addition, it is highly radiopaque and is recommended for all root canal preparation techniques.

→ optimal viscosity

Its optimal viscosity allows fast and simple use. Working time is 4 hours, allowing time intensive root canal treatment and corrections.

→ biocompatible

2Seal is biocompatible and has minimal shrinkage. Because of its solubility rate, **2Seal** is not affected by contact with bodily fluids over the years. It remains dense and does not allow the penetration of bacteria in the root canal.

Distributed in the U.S.A. by:



608 Rolling Hills Drive, Johnson City, TN 37604
800-992-7767

Ref. 1020 **2Seal**
Root canal sealing material

Contents:
4 ml paste A
4 ml paste B



Made in Germany

Rx Only

Manufactured for:

VDW GmbH
P.O.Box 830954
D-81709 Munich
Phone +49 89 62734-0 • Fax +49 89 62734-195
info@vdw-dental.com • www.vdw-dental.com



R-LIT2SEAL 8/05B

Endo Easy Efficient