

# Steel Vented Posts (SVP)



## Stainless Steel and Titanium Vented Post Kits

- Stainless steel posts are surgical grade for greater strength
- Parallel sides
- Non-engaging
- Comparable to Whaledent Para Posts™



Kits include 3 posts of 7 different sizes,  
1 each of corresponding drills and silicone stops

Size	mm	Steel Vented Posts	Drills	Titanium Posts	Drill Band Color
3	.90	SVP3STD	SVP3D3	TVP3STD	orange
4	1.00	SVP4STD	SVP4D4	TVP4STD	pink
4.5	1.12	SVP45STD	SVP45D45	TVP45STD	gray
5	1.25	SVP5STD	SVP5D5	TVP5STD	tan
5.5	1.37	SVP55STD	SVP55D55	TVP55STD	purple
6	1.50	SVP6STD	SVP6D6	TVP6STD	white
7	1.75	SVP7STD	SVP7D7	TVP7STD	blue

Posts available in refill  
packages of 10

Drills available in refill  
packages of 2

YOUR ROYDENT DISTRIBUTOR:



Johnson City, TN 37604  
800-992-7767 | www.roydent.com

R-VP 2/04

