

Gates Glidden Drills and Peeso Reamers



Gates Glidden Drills

- Enlarges coronal and middle portions of canal through bulk removal of tooth structure
- Sizes 1-6, identified by rings on shank
- Stainless steel instrument
- Flexible shaft
- 28mm (short) and 32mm (standard) lengths

Packaged 6 per pack

Peeso Size	32mm Standard	28mm Short
1	9991	999S1
2	9992	999S2
3	9993	999S3
4	9994	999S4
5	9995	999S5
6	9996	999S6
1-6	99916	999S16

Gates Glidden Size	32mm Standard	28mm Short
1	8881	888S1
2	8882	888S2
3	8883	888S3
4	8884	888S4
5	8885	888S5
6	8886	888S6
1-6	88816	888S16

Peeso Reamers

- For bulk removal of debris and gutta percha
- Creates post space
- Sizes 1-6, identified by rings on shank
- Stainless steel instrument
- Flexible shaft
- 28mm (short) and 32mm (standard) lengths

Packaged 6 per pack

