

Products for a Fast and Reliable Fill



A Look Inside:

- 2Seal and 2Seal easymiX
- Paste Carrier
- Gutta Percha
- Finger Pluggers
- Paper Points
- Finger Spreaders

2Seal and 2Seal easymiX



Take the Mess out of Mixing *Convenient auto-mix syringe with 360° Flex Tip*

- ✓ Automatic 1:1 mix ratio
- ✓ Mix tip creates less waste than manual mixing
- ✓ Dispenses only the small amount of sealer needed
- ✓ Consistently delivers controlled, precise mix
- ✓ Easy to use with slim adjustable 360° flex tip.
- ✓ Works with all obturation techniques; warm, cold and thermal.



The quality of your sealer is critical for a biocompatible and long-term tight root canal filling. Epoxy amine-based sealers have been thoroughly tested in more than 90 scientific studies and well-proven in clinical applications.

Excellent biocompatibility:

- Highly biocompatible
- Extruded material is easily resorbed and encapsulated by the surrounding tissue
- Antimicrobial


For all obturation methods:

- Recommended for warm and cold obturation techniques
- Optimal viscosity, reaches all canal areas
- Excellent radiopacit

Dense and durable:

- Outstanding flow characteristics
- Low shrinkage and low solubility
- High adhesion to gutta percha
- Long-term dimensional stability

360° Flex Tip allows for easier, more precise application.



2Seal and 2Seal easymiX



2Seal easymiX Starter Kit



2Seal easymiX Mixing Syringes



2Seal easymiX Mixing Tips

2Seal is available in a standard 2-tube pack and also as 2Seal easymiX in a convenient, double-chambered mixing syringe for easy and fast application. Available for order as kits, syringes only or mixing tips only for your convenience.

Product Specifications:

- Working time: 4 hours
- Setting time: 8 hours
- Radiopacity: 13.6 mm/mm Al
- Flow: 36 mm
- Film thickness: 26 µm
- Shrinkage: 1.76%
- Solubility: 0.31%

Item Number	Product	Description
102500R	2Seal easymiX Starter Kit	2Seal easymiX Mixing Syringe, 15 g 20 x 360° Flex Tips 1 Mixing Pad 1 Illustrated Tip Card 1 Organizer
102501R	Mixing Syringe Refill	2 x 2Seal easymiX Mixing Syringe, 15 g
102502R	Mixing Tips	40 x 360° Flex Tips




2Seal Tubes

2Seal is also available for order in tube packs which include one tube each of Paste A and Paste B plus a mixing pad tip.

Item Number	Product	Description
1020USR	2Seal – Tubes	1 x 4 ml Paste A 1 x 4 ml Paste B 1 x Mixing Pad Tips

Finger Plugger • Finger Spreader • Paste Carrier




Finger Plugger

- ✓ Pre-applied silicone stops, color coded by length
- ✓ Sizes correspond to standard K-Files and Reamers
- ✓ Plastic handles are color-code for easy use
- ✓ Superb tactile feel
- ✓ Highly polished to help prevent sticking
- ✓ Flat tip to aid vertical condensation

packaged 6 to a pack

Size	21mm [white stop]	25mm [yellow stop]
15/White	992115	992515
20/Yellow	992120	992520
25/Red	992125	992525
30/Blue	992130	992530
35/Green	992135	992535
40/Black	992140	992540
15-40	9921160	9925160



Finger Spreader

- ✓ Pre-applied silicone stops, color-coded by length
- ✓ Corresponds to standard K-Files and Reamers
- ✓ Handles are color-coded for easy use
- ✓ Superb tactile feel
- ✓ Highly polished to help prevent sticking
- ✓ Pointed tip to aid in lateral condensation

packaged 6 to a pack

Size	25mm [yellow stop]
15/White	952515
20/Yellow	952520
25/Red	952525
30/Blue	952530
35/Green	952535
40/Black	952540
15-40	9525160




Paste Carrier

- ✓ Manufactured from stainless spring steel
- ✓ Tightly wound safety spring at the base of the shaft
- ✓ Cone shaped to help air bubbles escape
- ✓ Color-coded shank

packaged 4 to a pack

Size	21mm	25mm
25/Red	932125	932525
30/Blue	932130	932530
35/Green	932135	932535
40/Black	932140	932540
25-40	9321190	9325190



Paste Carrier

Absorbent Paper Points



Absorbent Paper Points

.02 Taper ISO and Accessory Sized

- ✓ Uniformly machine-rolled
- ✓ Blunt tips for rapid absorption
- ✓ Super absorbent with no starches or fillers to leave residue in canals
- ✓ Ridged for easy insertion and removal
- ✓ Matched precisely to the dimensions of standardized instruments
- ✓ ISO color-coding to distinguish sizes easily
- ✓ Available in bulk or sterile packaging

.04 & .06 Taper ISO Sized

- ✓ Hand-rolled with color-coded ends
- ✓ Available in sterile cell packs
- ✓ Matches .04 & .06 rotary instruments

ISO Color-Coded

.02 Taper ISO (Cell, Sterile)

Size	29mm
15	0560015Z
20	0560020Z
25	0560025Z
30	0560030Z
35	0560035Z
40	0560040Z
45	0560045Z
50	0560050Z
55	0560055Z
60	0560060Z
70	0560070Z
80	0560080Z
15-40	0560160Z
45-80	0560165Z

Accessory Sized

(Cell, Sterile)

Size	28mm
XXFine	0565XXFZ
XFine	0565XFZ
Fine	0565FZ
Medium	0565MZ
Coarse	0565CZ
XCoarse	0565XCZ
XXFine-Coarse (Assorted)	0565XXFCZ

Accessory Sized (Bulk)

Size	28mm
XXFine	0555XXFZ
XFine	0555XFZ
Fine	0555FZ
Medium	0555MZ
Coarse	0555CZ
XCoarse	0555XCZ
XXFine-Coarse (Assorted)	0555XXFCZ

ISO Color-Coded Standardized

.02 Taper ISO (Bulk)

Size	29mm
15	0550015Z
20	0550020Z
25	0550025Z
30	0550030Z
35	0550035Z
40	0550040Z
45	0550045Z
50	0550050Z
55	0550055Z
60	0550060Z
70	0550070Z
80	0550080Z
15-40	0550160Z
45-80	0550165Z

ISO Color-Coded

.04 Taper ISO (Cell, Sterile)

Size	28mm
15	604015R
20	604020R
25	604025R
30	604030R
35	604035R
40	604040R
15-40	604160R

ISO Color-Coded

.06 Taper ISO (Cell, Sterile)

Size	28mm
15	606015R
20	606020R
25	606025R
30	606030R
35	606035R
40	606040R
15-40	606160R

Gutta Percha Points

Gutta Percha Points

.02 Taper ISO and Accessory Sized

- ✓ Excellent X-ray visibility
- ✓ Only pure, high-quality raw materials
- ✓ Uniformly machine-rolled with flattened ends
- ✓ Used in lateral, vertical or thermal condensation
- ✓ ISO color-coding to distinguish sizes easily
- ✓ Stiffness avoids bending in canals
- ✓ Ideal elasticity to adapt to canal walls
- ✓ Latex free

Increased Taper .04 and .06

- ✓ Hand-rolled with color-coded ends
- ✓ Allows for single cone obturation
- ✓ Excellent X-ray visibility
- ✓ Cadmium and latex-free
- ✓ Matches .04 & .06 rotary instruments



Gutta Percha Points .02 Taper ISO Color-Coded

Size	29mm
15	0525015Z
20	0525020Z
25	0525025Z
30	0525030Z
35	0525035Z
40	0525040Z
45	0525045Z
50	0525050Z
55	0525055Z
60	0525060Z
70	0525070Z
80	0525080Z
90	0525090Z
100	0525100Z
110	0525110Z
120	0525120Z
130	0525130Z
140	0525140Z
15-40	0525160Z
45-80	0525165Z
90-140	0525175Z

Gutta Percha Points Accessory Sized

Size	29mm
XFine	0530XFZ
XFine-Med(Asst'd)	0530FMZ
Fine	0530FZ
Fine-Fine	0530FFZ
Fine-Med	0530FMZ
Medium-Fine	0530MFZ
Medium	0530MZ
Medium-Large	0530MLZ
Large	0530LZ

Gutta Percha Points .04 Taper ISO Color-Coded

Size	29mm
15	504015R
20	504020R
25	504025R
30	504030R
35	504035R
40	504040R
15-40	504160R

Gutta Percha Points .06 Taper ISO Color-Coded

Size	29mm
15	506015R
20	506020R
25	506025R
30	506030R
35	506035R
40	506040R
15-40	506160R

Other available Roydent procedural catalogs include:

- Access
- Cleaning and Shaping
- Restoration
- Specialty
- Bite Registration & Impression Material

To order Roydent products, call 1.800.992.7767



608 Rolling Hills Drive
Johnson City, TN 37604
roydent.com • info@roydent.com
1.800.992.7767



OBT559/13