

discover

Aha!™

Make a Great First Impression—
Every Time.

Introducing the Aha! Family of Impression Materials

We understand how frustrating it can be to THINK you have the perfect impression only to have molds chip, break, or tear in processing. With Roydent's new Aha! family of impression materials, you can feel confident knowing you will get a solid, durable impression that will hold up and stay true. You get a great mold AND a happy patient, the FIRST time.

Now that truly IS a great impression!



SPECIAL INTRODUCTORY OFFER: Buy 2, Get 1 FREE!

Includes all Aha! products shown above | code RW1

ROYDENT™
DENTAL PRODUCTS

The boutique brand you've come to love & trust.

roydent.com | 1.800.992.7767 | [f](#) [y](#) [in](#)

*Redemption: To receive your FREE product, Fax your qualified dealer invoice dated 4/1/17 to 9/30/17 to: Roydent Dental Products (888) 769-3368 with Promo Code RW1 identified on or before 10/31/17. One qualifying invoice per redemption. Please indicate name, size and length of free product to be of equal or lesser value. Allow 4-6 weeks for delivery. Offers available to all Roydent customers. Not to be combined with any other Roydent offer. USA only.

Consistent, reliable impressions. Every time.



Aha!™ Fast Bite Registration Material

NEW SKU: K11961R

- Shore-hardness D 43
- Working time: 15 sec, Intraoral set time: 45 sec, Total set time: 1 min
- For “small applications” due to shorter working time
- Addition-curing elastomeric bite registration materials



Aha!™ Bite Registration Material

NEW SKU: K11932R

Bulk Package: K11932BR

- Shore-hardness D 43
- Working time: 30 sec, Intraoral set time: 90 sec, Total set time: 2 min
- For applications requiring a hard bite registration silicone
- Elastomeric bite registration materials intended for making accurate occlusal records



Aha!™ VPS Impression Material, X-Light

NEW SKU: K11892R

Bulk Package: K11892BR

- Working time: 2 min, Working time at 35°C/ 95°F: 1 min, Intraoral set time: 2 min, Total set time: 4 min
- Fabricating full or partial dentures
- Reline impressions
- Addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation



Aha!™ Putty Fast Set

NEW SKU: K11121R

- Working time: 2 min, Intraoral set time: 2 min, Total set time: 4 min
- Can be used for two-step or one-step putty-wash impression technique
- Can be used as two-step putty-wash impression technique using a foil (plastic putty spacer)
- Can be used as one-step putty impression technique for forming functional peripheries



Aha!™ VPS Heavy Body Impression Material | **NEW!**

NEW SKU: K13541R

- Working time: 1 min 20 sec, Intraoral set time: 2 min, Total set time 3 min, 20 sec
- Reduced final hardened
- Displays high initial hydrophilicity
- Addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation.

ROYDENT™
D E N T A L P R O D U C T S

The boutique brand you've come to love & trust.

roydent.com | 1.800.992.7767 |   