

K-Files

- Pre-applied silicone stops, color-coded by length: 19mm – blue, 21mm – white, 25mm – yellow, 28mm – red, 31mm – black
- Square cutting blades in smaller sizes and triangular cutting blades in larger sizes provide safety, flexibility and sharpness
- Contoured ribbed handle provides exceptional tactile control for reduced finger fatigue
- Conical, non-cutting tip follows canal curve, providing reduced risk of ledging
- Handles are numbered and color-coded with ISO symbol for easier identification

Packaged 6 to a box



Hedstrom Files

- Pre-applied silicone stops, color-coded by length: 21mm – white, 25 mm – yellow, 28mm – red, 31mm – black
- Outstanding flexibility and high torsion strength due to optimal geometry and design
- Custom designed Helix Angle for each size instrument
- Contoured ribbed handle provides exceptional tactile control for reduced finger fatigue
- Handles are numbered and color-coded with ISO symbol for easier identification

Packaged 6 to a box



K-Reamers

- Pre-applied silicone stops, color-coded by length: 21mm – white, 25mm – yellow, 28mm – red, 31mm – black
- Sharper angle positioning of the triangular profile produces a sharper blade
- Contoured handle provides exceptional tactile control for reduced finger fatigue
- Conical, non-cutting tip follows canal curve, providing reduced risk of ledging
- Handles are numbered and color-coded with ISO symbol for easier identification

Packaged 6 to a box





K-Files



Hedstrom Files



K-Reamers

Hedstrom Files

Size	21mm	25mm	28mm	31mm
08	07321008R	07325008R	07328008R	07331008R
10	07321010R	07325010R	07328010R	07331010R
15	07321015R	07325015R	07328015R	07331015R
20	07321020R	07325020R	07328020R	07331020R
25	07321025R	07325025R	07328025R	07331025R
30	07321030R	07325030R	07328030R	07331030R
35	07321035R	07325035R	07328035R	07331035R
40	07321040R	07325040R	07328040R	07331040R
08-10	07321155R	07325155R	-	07331155R
15-40	07321160R	07325160R	07328160R	07331160R
45	07321045R	07325045R	07328045R	07331045R
50	07321050R	07325050R	07328050R	07331050R
55	07321055R	07325055R	07328055R	07331055R
60	07321060R	07325060R	07328060R	07331060R
70	07321070R	07325070R	07328070R	07331070R
80	07321080R	07325080R	07328080R	07331080R
90	07321090R	07325090R	07328090R	07331090R
100	07321100R	07325100R	07328100R	07331100R
110	07321110R	07325110R	07328110R	07331110R
120	-	07325120R	-	07331120R
130	-	07325130R	-	07331130R
140	-	07325140R	-	07331140R
45-80	07321165R	07325165R	07328165R	07331165R
90-110	-	07325170R	-	07331170R

K-Files

Size	19 mm	21mm	25mm	28mm	31mm
06	06319006R	06321006R	06325006R	06328006R	06331006R
08	06319008R	06321008R	06325008R	06328008R	06331008R
10	06319010R	06321010R	06325010R	06328010R	06331010R
15	06319015R	06321015R	06325015R	06328015R	06331015R
20	06319020R	06321020R	06325020R	06328020R	06331020R
25	06319025R	06321025R	06325025R	06328025R	06331025R
30	06319030R	06321030R	06325030R	06328030R	06331030R
35	06319035R	06321035R	06325035R	06328035R	06331035R
40	06319040R	06321040R	06325040R	06328040R	06331040R
06-10	-	06321150R	06325150R	06328150R	06331150R
08-10	-	06321155R	06325155R	-	06331155R
15-40	06319160R	06321160R	06325160R	06328160R	06331160R
45	-	06321045R	06325045R	06328045R	06331045R
50	-	06321050R	06325050R	06328050R	06331050R
55	-	06321055R	06325055R	06328055R	06331055R
60	-	06321060R	06325060R	06328060R	06331060R
70	-	06321070R	06325070R	06328070R	06331070R
80	-	06321080R	06325080R	06328080R	06331080R
90	-	06321090R	06325090R	06328090R	06331090R
100	-	06321100R	06325100R	06328100R	06331100R
110	-	06321110R	06325110R	06328110R	06331110R
120	-	-	06325120R	-	06331120R
130	-	-	06325130R	-	06331130R
140	-	-	06325140R	-	06331140R
45-80	-	06321165R	06325165R	06328165R	06331165R
90-110	-	06321170R	06325170R	06328170R	06331170R
90-140	-	-	06325175R	-	06331175R

K-Reamers

Size	21mm	25mm	28mm	31mm
06	05321006R	05325006R	05328006R	05331006R
08	05321008R	05325008R	05328008R	05331008R
10	05321010R	05325010R	05328010R	05331010R
15	05321015R	05325015R	05328015R	05331015R
20	05321020R	05325020R	05328020R	05331020R
25	05321025R	05325025R	05328025R	05331025R
30	05321030R	05325030R	05328030R	05331030R
35	05321035R	05325035R	05328035R	05331035R
40	05321040R	05325040R	05328040R	05331040R
08-10	05321155R	05325155R	-	05331155R
15-40	05321160R	05325160R	05328160R	05331160R
45	05321045R	05325045R	05328045R	05331045R
50	05321050R	05325050R	05328050R	05331050R
55	05321055R	05325055R	05328055R	05331055R
60	05321060R	05325060R	05328060R	05331060R
70	05321070R	05325070R	05328070R	05331070R
80	05321080R	05325080R	05328080R	05331080R
90	05321090R	05325090R	05328090R	05331090R
100	05321100R	05325100R	05328100R	05331100R
110	05321110R	05325110R	05328110R	05331110R
120	-	05325120R	-	05331120R
130	-	05325130R	-	05331130R
140	-	05325140R	-	05331140R
45-80	05321165R	05325165R	05328165R	05331165R
90-110	05321170R	05325170R	05328170R	05331170R

YOUR ROYDENT DISTRIBUTOR:



Johnson City, TN 37604
800-992-7767 | www.roydent.com

R-ZipFiles 2/04

