

# MAP SYSTEM

MICRO - APICAL PLACEMENT





Introducing the  
Hooked Needle.  
Now Available In  
Four Sizes.



## Features of the MAP System

- Ideal for placement of endodontic repair cements such as ProRoot® MTA root canal repair material
- Available in Complete Endo and Intro Endo Kits
- Triple-angled needles for apical placement of materials
- Curved needles for coronal material placement
- Hooked needles for retrograde obturation, in sizes .9mm, 1.1mm, 1.3mm & 1.5mm
- Color Coding for easier tip identification
- Aluminum sterilization case
- Cleaning Kit included

**FOR RETROGRADE OBTURATION**  
after apical resection on anterior teeth**FOR RETROGRADE OBTURATION**  
after apical resection on both anterior and posterior teeth**FOR ORTHOGRADE OBTURATION**  
perforation, pulp capping, root-end filling**Item #****Description**

9105	Complete Endo Kit - 4 Triple-Angled Needles with 6 Red, 6 Yellow PEEK Pistons, 2 Curved Needles with 2 NiTi Pistons, Syringe, 2 Stainless Wells, 2 Cleaning Devices, 1 Set of Rings, 1 Ring Remover, 1 Sterilizable Alum. Case
9200	Intro Endo Kit - 2 Curved Needles with 2 NiTi Pistons, 2 Hooked Needles with 3 Red, 3 Yellow PEEK Pistons, Syringe, 1 Stainless Well, 2 Cleaning Devices, 1 Sterilizable Alum. Cases

**Refill Items**

9110	Stainless Steel Syringe (1)
9120	Triple-angled Needle - .9mm/right (1) - Yellow
9130	Triple-angled Needle - .9mm/left (1) - Yellow
9140	Triple-angled Needle - 1.1mm/right(1) - Red
9150	Triple-angled Needle - 1.1mm/left (1) - Red
9160	Curved Needle - .9mm (1) - Yellow
9170	Curved Needle - 1.1mm (1) - Red
9172	Hook Needle - .9mm (1) - Yellow
9174	Hook Needle - 1.1mm (1) - Red
9176	Hook Needle - 1.3mm (1) - Blue
9178	Hook Needle - 1.5mm (1) - Green
9180	Stainless Steel Mixing Well (1)
9190	Aluminum Sterilization Case - Large (1)
9195	Aluminum Sterilization Case - Small (1)
9205	Straight Needle - .9mm (1) - Yellow
9210	Straight Needle - 1.1mm (1) - Red
9215	Needle Cleaning Device - .9mm (1) - Yellow
9220	Needle Cleaning Device - 1.1mm (1) - Red
9222	Needle Cleaning Device - 1.3mm (1) - Blue
9224	Needle Cleaning Device - 1.5mm (1) - Green
9225	Peek Piston -.9mm (6) - Yellow
9230	Peek Piston - 1.1mm (6) - Red
9232	Peek Piston - 1.3mm (6) - Blue
9234	Peek Piston - 1.5mm (6) - Green
9235	NiTi Piston - .9mm (1) - Yellow
9240	NiTi Piston - 1.1mm (1) - Red
9245	Ring Remover (1)
9250	Silicone Rings - .9mm (50) - Yellow
9255	Silicone Rings - 1.1mm (50) Red

YOUR ROYDENT DISTRIBUTOR:

800.992.7767  
www.roydent.com

R-MAP 10/04